

CLINICAL EVIDENCE

VEXODERM DRESSING



**SKIN AND WOUND CARE PRODUCTS
WITH ANTIOXIDANT AND REGENERATIVE
PROPERTIES**

WWW.VEXODERM.COM

LACERATED-BRUISED WOUND

Colony cat, female. Adult.

Lacerated-bruised wound from rubbing.
Freshly contaminated wound.

Previous treatment

None.

Dr. Alessandro Negro, DVM, MSc

Soft tissue surgery and wound care.
Medical and surgical oncology.

Vexoderm treatment

Application of **Vexoderm Dressing**.

Week 1: cleaning of the wound bed and replacement of the dressing. Very lively granulation. Margins begin to contract.

Week 2: very evident wound contraction.

The volunteer proposes to apply **Vexoderm Gel+** every 8 hours to avoid further sedation.

Week 3: Tissues are closed to healing.

Without sedation, getting too close is risky.

Week 4: the wound is completely healed.



BEGINNING



WEEK 1



WEEK 2



WEEK 3



WEEK 4

EXPLODED ABSCESS

**European Shorthair cat, sterilised female.
9 years old.**

Exploded abscess in the sacro-ischial region.

Previous treatment

Attempt to treat it with sugar.
Sedation, trichotomy, washing and cleansing.
Disinfection with chlorhexidine and povidone iodine.
Surgical debridement and buttonholes for tie overs.

Vexoderm Dressing treatment

Application of Vexoderm Dressing and tie-over. Excellent granulation and contraction of margins.

Dr. Alessandro Negro, DVM, MSc

Soft tissue surgery and wound care.
Medical and surgical oncology.



BEGINNING



WEEK 1



WEEK 2

TRAUMATISM

Dog, male.

Exposed bone due to traumatism.

Previous treatment

None.

Vexoderm Dressing treatment

No. of changes: 4 (once a week).

No. of dressings: 4 Vexoderm Dressing 9x4.

VetCenter Clinic, Bucharest.

Dr. Isabelle Caloian.



WEEK 1



WEEK 2



WEEK 3



WEEK 4

TRAUMATIC WOUND

Dog, male. Fox Terrier. 10 years old.

Traumatic wound in one leg.

Previous treatment

None.

Vexoderm Dressing treatment

No. of changes: 3 (once a week).

No. of dressings: 3 Vexoderm 9x4.

Wound closure: week 3.



BEGINNING



WEEK 1



WEEK 2



WEEK 3

TRAFFIC ACCIDENT

Dog, female. Border Collie. 11 months old.

Traumatic wound in one leg.

Previous treatment

None.

Vexoderm Dressing treatment

No. of changes: 3 (once a week).

No. of dressings: 3 Vexoderm 9x4.

Wound closure: week 3.



BEGINNING



WEEK 1



WEEK 2



WEEK 4



WEEK 5



WEEK 6

TUMOUR RESECTION

Dog, male. Golden Retriever. 9 years old.
Surgical resection of hepatoid gland tumour in the tail.

Previous treatment
Antibiotics.

Vexoderm Dressing treatment

No. of changes: 2 (once a week).
No. of dressings: 2 Vexoderm 9x4.
Wound closure: week 2.

Prevents potential wound dehiscence.
No maceration.



BEGINNING



WEEK 1



WEEK 2



WEEK 11

NECROTIC ULCER

**Dog, female. American Staffordshire.
2 years old.**

Suppurative necrotic ulcer. Infected.

Previous treatment

Antibiotics / NSAIDs. Surgical debridement.

Vexoderm Dressing treatment

No. of changes: 7 (first change after 3 days and then once a week).

No. of dressings: 7 Vexoderm 9x4.

Wound closure: week 7.



BEGINNING



BEGINNING



WEEK 1



WEEK 2



WEEK 3



WEEK 4



WEEK 7



WEEK 40

BURN

Cat, female. Felis Vulgaris. 10 months old.

Burn caused by deliverately putting the cat in the oven.

Previous treatment



Moisturising products. Antibiotics. Surgical debridement.

Vexoderm Dressing treatment

No. of changes: 9 (two changes the first week and then once a week).

No. of dressings: 12 Vexoderm 10x7.
2 Vexoderm 9x4.


Wound closure: week 8.



BEGINNING

WEEK 1

WEEK 2



WEEK 3

WEEK 4

WEEK 6

WEEK 6

TRAUMA

Cat, female. European Shorthair cat. 1 year old.

Tissue loss due to trauma.
10 days in evolution.

Previous treatment

Antibiotics. Saline cures.

Vexoderm Dressing treatment

No. of changes: 5 (once a week).

No. of dressings: 5 Vexoderm 9x4.

Wound closure: week 5.



BEGINNING



WEEK 1



WEEK 2



WEEK 3



WEEK 4



WEEK 5



WEEK 6

POLYTRAUMATISM

Cat, male. 7 years old.

Polytraumatism, necrotic wound, emphysema and edema.

Previous treatment

Suture. Spray (silver).

Vexoderm Dressing treatment

No. of changes: 10 (once a week).

No. of dressings: 10 Vexoderm 9x4.

Wound closure: week 10.



BEGINNING



WEEK 1



WEEK 2



WEEK 3



WEEK 5



WEEK 9



WEEK 9



WEEK 10

CUT WITH HYPERGRANULATION

Horse, male. Crossbreed. 15 years old.

Cut, hypergranulation "Proud Flesh".

Previous treatment

Antibiotics and cures. Surgical debridement.

Vexoderm Dressing treatment

No. of changes: 32 (once a week).

No. of dressings: 24 Vexoderm 10x7
and 16 Vexoderm 9x4.

Wound closure: week 32.

Tissue resection at weeks -1, 2, 4 & 5.



BEGINNING



WEEK 1



WEEK 4




WEEK 10



WEEK 17



WEEK 28



**6 MONTHS
AFTER
CLOSURE**

CHRONIC WOUND

Horse, male. Thoroughbred.
2 years old.

Chronic wound. 5 months in evolution.

Previous treatment

Antibiotics 2 days before Vexoderm.

Vexoderm Dressing treatment

No. of changes: 5 (every 7-10 days).

No. of dressings: 5 Vexoderm 10x7.

Wound closure: week 7.



BEGINNING



WEEK 2



WEEK 3



WEEK 4



WEEK 6



WEEK 7



WEEK 12

CHRONIC ULCER

Horse.

Chronic ulcer of unknown ethiology.

Previous treatment

Multiple different treatments including honey.

Vexoderm Dressing treatment

No. of changes: 7 (once a week).

No. of dressings: 4 Vexoderm 10x7.

6 Vexoderm 9x4.

Wound closure: treatment with Vexoderm ended before wound closure (week 7).



BEGINNING



WEEK 1



WEEK 2



WEEK 3



WEEK 4



WEEK 6



WEEK 7

CUT WITH FLAP

**Horse, female. Thoroughbred.
22 years old.**

Cut with flap. 5 months in evolution.

Previous treatment

Antibiotics and anti-inflammatory drugs.

Vexoderm Dressing treatment

No. of changes: 7 (once a week).

No. of dressings: 7 Vexoderm 10x7.

Wound closure: week 7.



BEGINNING



FLAP REMOVAL



WEEK 2



WEEK 3



WEEK 5



WEEK 7

TRAUMA

Dog, female. 3 months old.
Tissue loss due to traumatism.


Previous treatment

Manuka honey plus antimicrobial cream.

Vexoderm Dressing treatment

Vexoderm Dressing 10x7, 9x4
and Vexoderm Gel+.

Vet4Life Clinic, Bucharest.



BEGINNING



DAY 7



DAY 7



DAY 13

ULCER

Dog, female. Bulldog. 4 years old.

Ulcer due to laparotomy.
Possible reaction to suture.

Previous treatment

Systemic antibiotics.
Dressing 6 days after surgery.

Vexoderm Dressing treatment

Wound activation and healing.
Treatment with Vexoderm Dressing
until wound closure.



BEGINNING



WEEK 2



WEEK 3

CHRONIC SORE

Quarter Horse, male. 12 years old.



Chronic saddle sore in both sides for over 5 months.

Previous treatment
Ointment.

Vexoderm Dressing treatment

Vexoderm Dressing is applied with dressing changes once a week for two weeks until complete closure of the wounds.

A few weeks later hair growth is observed.

**BEGINNING****WEEK 1****WEEK 11**